


CAPILLARYS **3**
DBS

*The perfect match for screening
hemoglobin disorders
in newborns*



sebia 

The new language of life

High throughput solution for newborn hemoglobin screening



Newborn Hemoglobin screening has never been so easy...

SAMPLING

PUNCHER

ELUTION

ANALYSIS

RESULT

Main benefits

1
Ease of use, speed and reduced hands-on time

2
Automatic result flags let you **focus on pathological samples**

3
Clear detection of **Hemoglobin variants** (S, C, D, E...) as well as **Thalassemia (Hb Bart's)**

Key features



High throughput instrument

70 samples/hr



Full traceability

From sample punch to result



High autonomy

8 full 96-well microplates could be analyzed in a single run



Automatic Reagent Control

RFID for all reagents and Hb AF control



Excellent Resolution

High voltage & temperature-controlled migration




Automatic presumptive identification of Hemoglobin Pattern


Quick check of pathological samples versus non-pathological samples (color-coded) due to automatic software interpretation

Results examples


FA + Bart's



FS (Sickle Cell Anemia)



FSC (Sickle Cell Disease)



4

Early detection of hemoglobin disorders allows appropriate treatment and follow up, **improving baby's life** while reassuring parents



Ordering information

CAPILLARYS 3 DBS

Instrument

• CAPILLARYS 3 DBS1244

Reagents

• CAPI 3 NEONAT Hb2506

• CAPI 3 NEONAT Hemolysing solution2563

Controls/ Calibrators

• CAPI 3 NEONAT Hb AF Control4757

• Hb AFSC Control4792


Consumables

• CAPILLARYS 3 DBS Microplates
(10 units - 96 well barcoded microplates)2578

• CAPILLARYS 3
CAPICLEAN2060

• CAPILLARYS 3
Wash Solution2062

• CAPILLARYS 3
Reagent Cups (24x14 units)2582



DISTRIBUTED BY



Bojana Uzelac

Product Manager

+381 63 531 743

bojana.uzelac@east-diagnostics.com

Mike Alasa 19, 11000 Belgrade
Republic of Serbia

east-diagnostics.com