


# Large menu of gel electrophoresis *on HYDRASYS 2 Scan Focusing*

**Semi-Automated  
Multi-Assay Analyzer  
Proven Technology**




**sebia**   
The new language of life

# Main Strengths

HYDRASYS 2 <sup>Scan Focusing</sup> gel electrophoresis instrument provides your laboratory with high-quality results.

Combining the HYDRASYS 2 <sup>Scan Focusing</sup> with our front-end pipetter, ASSIST, improves workflow and traceability.



## High Quality Result

- Have confidence when reporting your patients' results with gel electrophoresis assays.

## Multi-Functional Analyzer

- Ability to migrate, stain, & scan gels on one system.



## Powerful Software Solution

- PHORESIS Core Software provides tools for interpretation, data & result management and networking.

## Workflow Traceability

- Automation and standardization of samples preparation with ASSIST.

## All-in-One System Gel Electrophoresis


- All-in-one system performing migration, staining & scanning
- Small footprint fits any size work area
- Intuitive touchscreen interface

## Gel Automation with the ASSIST


- Preparation and automatic dispensing of samples from primary tube and antisera for dynamic mask application
- Auto-dilutes patient samples for IFE and CSF testing
- Traceability of primary tubes with positive sample ID with PHORESIS




# Multi gel sizes to fit your Laboratory's needs



**Protein (Serum)**  
PROTEIN(E) and  $\beta 1 - \beta 2$   
Gels sizes - 7, 15, 30, 54





**Immunofixation (Serum/Urine)**  
IMMUNOFIXATION  
Gels sizes - 1, 2, 4, 9




**BENCE JONES**  
Gels sizes - 1, 2, 4




**Hemoglobin**  
HEMOGLOBIN(E) ACID(E) & Alkaline  
Gels sizes - 7, 15




**HYDRASHIFT**  
daratumumab/isatuximab  
Immunofixation  
Gels sizes - 2, 4




**URINE PROFIL(E)**  
Immunofixation  
Gels sizes - 2, 4




**Iso-Electro Focusing - CSF & A1AT**  
CSF ISOFOCUSING  
Gels sizes - 3, 9




**Alpha 1 - Antitrypsin**  
Gels sizes - 18




**Enzymes**  
Isoforms: Lactate Dehydrogenase (LD),  
Creatine Kinase (CK) & Phosphatase  
Isoenzymes (ISO-PAL)  
Gels sizes - 7, 15



**Lipids**  
LIPOPROTEIN(E) - LIPO + Lp(a)  
Gels sizes - 7, 15, 30



**von WILLEBRAND MULTIMERS**  
Immunofixation  
Gels sizes - 5, 11



**APO E\***  
ISOFOCUSING  
Gels size - 18

## LARGE RANGE OF ASSAYS

- Serum & Urine Protein Electrophoresis (SPE/UPE)
- Serum & Urine Immunofixation (IFE, BENCE JONES)
- HYDRASHIFT daratumumab/isatuximab
- CSF and Alpha-1 Antitrypsin by ISOFOCUSING

- von WILLEBRAND MULTIMERS
- Iso-Enzymes, B2 transferrin,
- LDL/HDL cholesterol, lipoproteins
- HEMOGLOBIN(E) Electrophoresis
- APO E\*

\*For research use only

# Ordering information

## Instruments

• HYDRASYS 2 Scan Focusing .....	1202
• ASSIST .....	1218

## Reagents

• HYDRAGEL 7/15/30 PROTEIN(E) .....	4100/4120/4140
• HYDRAGEL 54 PROTEIN .....	4145
• HYDRAGEL 7/15/30 PROTEIN b1b2 .....	4101/4121/4141
• HYDRAGEL 54 b1b2 .....	4146
• HYDRAGEL 7/15 HR .....	4102/4122
• HYDRAGEL 5 PROTEINURIE .....	4115
• HYDRAGEL 6/12 IF PENTA .....	4341/4342
• HYDRAGEL 1/2 IF .....	4801/4802
• HYDRAGEL 4/9 IF .....	4804/4809
• HYDRAGEL 2 IF/BJ HR .....	4803
• HYDRAGEL 1 BENICE JONES .....	4821
• HYDRAGEL 2/4 BENICE JONES .....	4822/4824
• HYDRAGEL 2/4 URINE PROFIL(E) .....	4831/4832

\* Non-exhaustive list, contact your local representative for more information about the Sebia gel range

\*\* For research use only

## Reagents

• HYDRASHIFT 2/4 daratumumab (20 tests) .....	4639
• HYDRASHIFT 2/4 isatuximab (40 tests) .....	4642
• HYDRAGEL 5/11 von WILLEBRAND MULTIMERS .....	4359/4360
• HYDRAGEL 3/9 CSF ISOFOCUSING .....	4353/4355
• HYDRAGEL 18 A1AT ISOFOCUSING .....	4357
• HYDRAGEL 15/30 LIPO + Lp(a) .....	4124/4144
• HYDRAGEL 15/30 LIPOPROTEIN(E) .....	4134/4138
• HYDRAGEL 15/30 LDL/cholesterol .....	4135/4139
• HYDRAGEL 7/15/30 ISO-LDH .....	4110/4130/4136
• HYDRAGEL 7/15/30 ISO-CK .....	4111/4131/4137
• HYDRAGEL 7/15 ISO-PAL .....	4112/4132
• HYDRAGEL 6 b2 TRANSFERRIN(E) .....	4853
• HYDRAGEL 18 APO E ISOFOCUSING ** .....	4349

DISTRIBUTED BY



**Aleksandra Mihajlov**

Product Manager

+381 63 531 743

[aleksandra.mihajlov@east-diagnostics.com](mailto:aleksandra.mihajlov@east-diagnostics.com)

Mike Alasa 19, 11000 Belgrade  
Republic of Serbia

[east-diagnostics.com](http://east-diagnostics.com)