

# **FastGene® PCR Adhesive PCR Foil**

- 138 x 79 mm (with edge), 118 x 79 mm (without edge)
- Suitable for Real-Time PCR applications
- Suitable for PE, PS and PP plates
- Prevents evaporation during PCR or storage
- Without sticky residue after the seal is peeled off
- End tabs for easy removal
- Resistant to DMSO
- Can be used at temperatures from -80 °C to +120 °C



Cat. No.: FG-93AC2 (100 sheets)

## Ordering information

Cat. No.	Product	Content
FG-93AC2	FastGene® Adhesive PCR Foil	100 sheets

# **FastGene® PCR Adhesive Seal Aluminium**

- 141 mm x 79 mm (with edge) and 121 mm x 79 mm (adhesive area)
- Soft aluminium for high conformability
- For use between -70 °C and +100 °C
- Adheres well to a wide range of plate material (PE, PS, PP)
- Easily piercable
- Pressure sensitive acrylate adhesive on the other side



Cat. No.: FG-93AC2 (100 sheets)

## Ordering information

Cat. No.	Product	Content
FG-93AF	FastGene® Adhesive Seal Aluminium	100 sheets

## Free Sample?

Just give us a call or write us an email if you would like to test our seals and receive a free sample!



+49 2421 554960



info@nippongenetics.de