



Pipet-Lite XLS+ Performance You Can Feel

Rainin Pipet-Lite™ XLS+ raises the bar on comfort, performance and control. Many volume sizes available in both LTS  and universal  formats.

Smooth operation

The sealing system on XLS+ pipettes delivers a plunger stroke that is smooth, precise and requires less force than before. Perfect for pipetting applications that require delicate control!

Highly accurate

Built to last and manufactured with premium quality materials, XLS+ pipettes are designed to deliver superior performance, year after year.

Autoclavable

XLS+ single-channel pipettes are designed for easy autoclaving where it counts: the complete liquid end, including the seals and piston. Autoclaving won't compromise performance or feel.

Total comfort

With sure-fit handles, contoured fingerhooks, lighter springs and low-drag seals, you can comfortably use Pipet-Lite XLS+ pipettes all day.

A better tip ejector

The tip ejector is made from a high density polymer that is corrosion free, autoclavable and highly resistant to chemical and physical insults. The ejector slides off in a single easy motion for cleaning.



All Rainin XLS+ pipettes are equipped with RFID tags for enhanced security, state-of-the-art inventory control and calibration tracking.

Secure volume lock

The snag-proof lock design on XLS+ pipettes prevents accidental volume drift while allowing for rapid volume adjustment, even while wearing gloves.

Advanced asset management

The RFID tag on board every Rainin XLS pipette stores calibration data that is electronically read by the Rainin SmartStand and EasyDirect software. See pages 57-58 for details.



GLP/GMP security

Minimize pipetting SOP deviations by locking protocol settings. Onboard service data is extremely secure and quick to access. Avoid calibration compliance errors by password protecting access to service interval alarms.



LTS	CATALOG No.	VOLUME RANGE
E4 XLS+	17014484	0.5–10 µL
	17014487	2–20 µL
	17014483	10–100 µL
	17014486	20–200 µL
	17014488	20–300 µL
	17014482	100–1000 µL
	17014485	200–2000 µL
	17012312	500–5000 µL
	17012313	1–10 mL
	17012314	2–20 mL

Rainin LTS pipettes are designed for use with LTS tips. Li-Ion battery and wall power supply are included.



UNV	CATALOG No.	VOLUME RANGE
E4 XLS+	17014491	0.5–10 µL
	17014494	2–20 µL
	17014490	10–100 µL
	17014493	20–200 µL
	17014495	20–300 µL
	17014489	100–1000 µL
	17014492	200–2000 µL
	17012353	500–5000 µL
	17012354	1–10 mL

Universal-Fit pipettes use universal tips. Li-Ion battery and wall power supply are included. For pipette specifications, see page 32.



Specialty Tips Unique Tips for Special Applications



Extended-length tips
For tall narrow tube work. BioClean Ultra.



Wide-orifice tips
Ideal for delicate samples. BioClean Ultra.



ShaftGuard™ 10µL tips
Avoid cross-contamination. BioClean Ultra.



Gel-Well™
Load gels in thin wells.

Pipet-Lite XLS + Multichannels

Precision, Performance and Comfort

The Rainin multichannel Pipet-Lite™ XLS™ + is one of the most technically and ergonomically advanced manual multichannel pipettes on the market.

Absolute consistency

New mechanical designs on all XLS+ multichannel pipettes assure the highest performance and industry-leading channel-to-channel consistency.



Fast tip loading

LTS tips for a perfect seal with minimal force, making them as easy to eject as they are to load. Absolutely no rocking or hand-tightening required!



UNV

CATALOG No.	VOLUME RANGE
1701 4413	0.1–2 µL
17014409	0.5–10 µL
17014412	2–20 µL
17014408	10–100 µL
17014411	20–200 µL
17014414	20–300 µL
17014407	100–1000 µL
17014410	200–2000 µL
17011801	500–5000 µL
17011795	1–10 mL

Pipet-Lite XLS+



All Rainin XLS+ pipettes are equipped with RFID tags for enhanced security, state-of-the-art inventory control and calibration tracking.

RAININ

Pipetting 360°



 **East Diagnostics**

Ovlašćeni distributer:

East Diagnostics d.o.o. Beograd
Mike Alasa 19, 11000 Beograd

Informacije:

milos.dzebric@east-diagnostics.com
+381 69 223 04 20

www.east-diagnostics.com

RAININ

Pipetting 360°



Distributed by

 **East Diagnostics**