


# Evolution makes the difference

CAPILLARYS  
OCTA **3**

Increasing your lab efficiency:  
CAPILLARYS 3 OCTA replaces CAPILLARYS 2,  
the most well-known capillary instrument worldwide



**sebia**

The new language of life

## Instrument

- Capillary electrophoresis instrument with 8 capillaries in parallel
- Instrument management via LCD touch screen
- Automatic or manual startup, shutdown & maintenance
- Automatic dilution segment management
- Automatic technique change using dedicated switch racks

## Assays & Throughput

- PROTEIN(E) 6 (serum) - 79 tests/h
- PROTEIN(E) 6 (urine) - 70 tests/h
- IMMUNOTYPING (serum or urine) - 9 tests/h
- Hb A1c - 43 tests/h
- HEMOGLOBIN(E) - 45 tests/h
- CDT - 49 tests/h

## Specifications

- Length: 90 cm
- Height: 54 cm
- Depth: 67 cm
- Weight: 75 kg

## Samples

- Serum, urine, whole blood
- Full traceability from the primary tube to final result
- Cap piercing and sample mixing for whole blood sampling on primary capped tubes for (Hb/Hb A1c)
- Initial simultaneous loading of 120 sample tubes

## Reagents

- Four flexible positions for reagents with automatic switch between vials
- Secondary reagent platform with cooled positions for antisera and additional reagents
- Full traceability with RFID tags on all reagent bottle and antisera

## Software

- PHORESIS
  - One network solution to connect all Sebia instruments around a single database
  - Expert software with specific features to manage your capillary menu

## Ordering information



DISTRIBUTED BY



**Aleksandra Bogdanović**

MSc product specialist  
+381 69 11 55 952

[aleksandra.bogdanovic@east-diagnostics.com](mailto:aleksandra.bogdanovic@east-diagnostics.com)

Mike Alasa 19, 11000 Belgrade  
Republic of Serbia

[east-diagnostics.com](http://east-diagnostics.com)